Strattice™ Reconstructive Tissue Matrix
Laparoscopic / Laparoscopique

**Instructions
For Use**

Surgical Mesh

Table of Contents

Instructions for Use English.....	3-6
Instructions for Use French.....	7-10

Device Description

Strattice™ Reconstructive Tissue Matrix (Strattice™ TM or Strattice™ surgical mesh) is a surgical mesh that is derived from porcine skin and is processed and preserved in a patented aqueous phosphate buffered solution containing matrix stabilizers. Strattice™ TM is intended to perform as a surgical mesh for soft tissue repair while presenting a scaffold to the patient's tissue. The structural properties minimize tissue adhesions to the mesh. The Strattice™ surgical mesh consists of a sterilized sheet of processed porcine dermis provided in prescribed different sizes, dimensions, and thicknesses and packaged in a double foil pouch.

Use of the surgical mesh provides for a strong and biocompatible implant and will incorporate into the patient's tissue with associated cell and microvascular ingrowth.

Animal studies show a low incidence in adhesion to the Strattice™ surgical mesh based on observation of minimal visceral tissue attachment.

Indications

Strattice™ surgical mesh is intended for use as a soft tissue patch to reinforce soft tissue where weakness exists and for the surgical repair of damaged or ruptured soft tissue membranes. Indications for use include the repair of hernias and/or body wall defects which require the use of reinforcing or bridging material to obtain the desired surgical outcome during open or laparoscopic procedures.

Strattice™ surgical mesh is supplied sterile and is intended for single patient one-time use.

Contraindications

- This surgical mesh is derived from a porcine source and should not be used in patients with known sensitivity to porcine material.
- Polysorbate 20 is a component of the phosphate buffered aqueous solution and therefore Strattice™ surgical mesh should not be used in patients with a known sensitivity to this material.

Warnings

- **Do not resterilize.**
- Do not use if the package is opened or damaged. Do not use if seal is broken or compromised.
- After use, handle and dispose of all unused product and packaging in accordance with accepted medical practices and applicable national and regional environmental laws on disposal of packaging and biological materials.
- Strattice™ surgical mesh cannot be reused once it has been removed from the packaging and/or is in contact with a patient without increased risk of patient-to-patient contamination and subsequent infection.

Precautions

- Discard surgical mesh if handling has caused possible damage or contamination.
- Discard surgical mesh if it is past the use-by-date of the product (indicated as 4 digit year, 2 digit month, and 2 digit day [YYYY-MM-DD]).
- Ensure that the surgical mesh is put into a sterile basin and covered with room temperature sterile saline or lactated Ringer's solution for a minimum of 2 minutes prior to implantation.
- Place surgical mesh in maximum possible contact with healthy, well-vascularized tissue to promote cell ingrowth and tissue remodeling.
- The Strattice™ surgical mesh should be hydrated and moist when the package is opened. If the Strattice™ surgical mesh is dry, do not use.
- In significantly contaminated or infected cases, utilize bioburden-reducing techniques to minimize contamination levels at the surgical site, including, but not limited to, appropriate drainage, debridement, negative pressure therapy, and/or antimicrobial therapy prior and in addition to implantation of the Strattice™ surgical mesh. (Presence of a significant microbial load may impact overall performance of the surgical mesh.)
- In large abdominal wall defect cases where midline fascial closure cannot be obtained, with or without separation of components techniques, utilization of Strattice™ surgical mesh in a bridged fashion is associated with a higher risk of hernia recurrence than when used to reinforce fascial closure.
- Refrain from using excessive force if inserting the mesh through a trocar.

Adverse Events

Adverse events are those typically associated with surgical mesh materials and their implantation procedures. Reported adverse events for Strattice™ surgical mesh used in hernia/body wall repair procedures include:

- Infection
- Hematoma
- Seroma formation
- Recurrence of tissue defect
- Fistula formation
- Inflammation
- Adhesion formation

Storage

- The Strattice™ surgical mesh is a sterile medical device that should be stored in a clean, dry location at -8°C to 30°C.
- It is to be stored in its original packaging.
- The use-by-date of the product is indicated as 4 digit year, 2 digit month, and 2 digit day (YYYY-MM-DD).

Sterilization

This product is sterilized by E-Beam.

Instructions For Preparing Strattice™ TM For Surgical Use

These instructions are designed to serve only as a general guideline. They are not intended to supersede institutional protocols or professional clinical judgment concerning patient care. Physicians should be familiar with surgical procedures and techniques involving biological surgical mesh before using Strattice™ surgical mesh.

Required Materials

- Sterile forceps
- Soaking fluid: room temperature sterile saline or room temperature sterile lactated Ringer's solution.
- One sterile basin per piece of the Strattice™ surgical mesh

Preparation

1. Open the carton and remove the foil pouch.
2. Peel open the outer foil pouch and remove the inner foil pouch using aseptic technique. The inner foil pouch is sterile and may be placed directly into the sterile field.
3. Open the inner pouch carefully and aseptically remove the surgical mesh. Always use sterile gloves or forceps when handling Strattice™ surgical mesh.
4. Place the surgical mesh in a sterile basin containing room temperature sterile saline or room temperature sterile lactated Ringer's solution for a minimum of 2 minutes.
5. Keep the surgical mesh in the room temperature sterile solution until ready for implantation. Surgical mesh can be retained in sterile solution for a maximum of 4 hours.

Implantation

1. Prepare the surgical site using standard techniques.
2. Strattice™ surgical mesh may be folded, trimmed or cut as required to fit the surgical site using aseptic technique, ensuring allowance for overlap.
3. Transfer Strattice™ surgical mesh to the surgical site using sterile gloves or forceps.
4. Suture and tack Strattice™ surgical mesh into place.

NOTE: Tension and suture placement are application dependent. For hernia repair applications, surgical experience with biological surgical mesh indicates that suturing the Strattice™ TM under physiologic tension with a minimum of 3cm–5cm overlap, or as much as required to reach healthy adjacent tissues, may produce improved outcomes. Use of permanent sutures is recommended. The choice of fixation device and fixation method should be determined by surgeon preference and the nature of the reconstruction to provide for adequate tissue fixation.

5. Complete the standard surgical procedure.
6. Discard any unused portions of Strattice™ surgical mesh as per institutional procedures.

Definitions

FIRM

DO NOT USE
IF PACKAGE IS
DAMAGED

OPEN HERE

CONTENTS

INNER POUCH

OUTER POUCH

CORRECT
SETTING

*INCORRECT
SETTING

Customer Support

*Call customer support on day of receipt for immediate return authorization.

For product complaints and potential adverse events, please contact your local Sales Representative or Customer Support.

This product and certain methods are covered by patents and patents pending, including: US 8,905,826; US 8,735,054; US 8,469,779; US 8,323,352; US 8,007,531; and US 7,476,249.

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© LifeCell Corporation 2015
LifeCell Corporation